Distributor:
Maxxspine Innovative GmbH
Große Hub 3a
65344 Eltville
Germany

D.A.R.T.

CEMENT INJECTION CANNULA FOR PEDICLE SCREWS
CANNULA INIEZIONE CEMENTO PER VITI PEDUNCOLARI
KANÜLE ZUR INJEKTION VON ZEMENT FÜR PEDIKELSCHRAUBEN
CÁNULA PARA LA INYECCIÓN DE CEMENTO PARA TORNILLOS PEDICULARES
CANULE INJECTION CIMENT POUR VIS PÉDICULAIRES

ENGLISH	CEMENT INJECTION CANNULA FOR PEDICLE SCREWS
---------	---

Instructions for use

- Before use check that the device package is not damaged. Check that the device is not damaged or kinked; if so, inform the manufacturer or your area representative. Take the device out of the package operating under aseptic conditions.
- Sedate and position the patient. The entire procedure must be performed in real-time under combined CAT and fluoroscopic guidance. Digital angiography may also prove useful. Perform local or general anaesthesia.
- Perform the necessary skin disinfection procedures. Prepare the various instruments contained in the kit required for vertebral fixation.
- Open the pedicles using a dedicated instrument and following the manufacturer's instructions for use.
- Insert the screws in the pedicles and position them following the procedures and instructions for use indicated by the manufacturer.
- Clean the canal with the dedicated wire contained in the kit (optional).
- Prepare and mix the cement according to the procedures indicated by the manufacturer.
- Fill the cannula with cement by connecting it to a dedicated syringe until it flows from the distal end.

- Insert this cannula in the vertebral pedicle and into the fixation screw you intend to fill.
- Let the cement flow from the cannula by exercising pressure from above with the dedicated pusher stylet.

Warnings

- Do not use the instrument if the package is damaged.
- Do not use the instrument beyond the expiry date.
- The device is single-use and must be destroyed after use. Reuse is strictly prohibited. Reuse or resterilization may compromise the structural integrity of the device and/or cause damage to the device which may cause injury, disease or death of the patient. Reuse or resterilization may also create a risk of device contamination and/or cause infection of the patient or cross-infection, including but not limited to transmission of infectious diseases from one patient to the next.
- Check that the screw is properly filled.

Precautions for use

The system may exclusively be used by physicians. After use, the device may constitute a risk of contamination and/or infection. Therefore, handle with care and dispose of the product in accordance with the regulations in force and medical practice. In order to minimise the risks or other complications associated with use of the cement, the following is recommended:

- The device should only be used by qualified physicians appropriately trained to perform the procedure.
- Use only cements intended for the fixation procedure and carefully read the instructions for use of the cement.
- Monitor the procedure using high-quality imaging devices in order to detect cement leakage or seepage from the application areas.
- It is also recommended not to use the device if particularly risky manoeuvres are necessary where it will be subjected to high torsional or bending forces which might compromise the integrity and proper functioning of the system.

Indications

The device is to be used for injecting cement into pedicle screws. It is a particularly useful instrument in the vertebral stabilization approach using cannulated and/or windowed pedicle screws. Thanks to this device, the PMMA is injected through the screw for better screw anchorage to the pedicle.

Pathological indications:

- Degenerative discopathy
- Fractures
- Tumours
- Osteoporosis

To be used under medical supervision only.

Packaging

The system is individually packed in laminated film and medical-grade paper bags and is sterilized by ethylene oxide (EtO) in compliance with ISO 11135 requirements. The device is sold in packs containing a number of pieces corresponding to that indicated on the label. The device is not reusable. If it is reused, there is a risk of altering device performance and/or the possibility of cross-contamination. Do not use the device if the package is open or damaged or contaminated.

Biopsybell declines all responsibility for any damage caused by improper use or different from what is indicated in this instruction leaflet.

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